EN

This medical device and Instructions for Use are intended for use by a lay person.

Read the instruction for use before starting the treatment with Sorbact® Compress.

Device Description

Sorbact® Compress is a DACC™-coated, non-absorbent bacteria and fungi binding wound contact layer, based on Sorbact® Technology, that allows passage of wound exudate into a secondary dressing. Sorbact® Compress prevents and treats wound infections and facilitates the wound healing process.

Intended Purpose

Sorbact® Compress is intended for use on exuding wounds such as surgical wounds and acute wounds.

Who can use Sorbact® Compress?

Sorbact® Compress is intended to be used on children and adults. The dressing can be used by healthcare professionals and lay persons in healthcare facilities and home environment.

How does Sorbact® Compress work?

Sorbact® Compress binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus* spp., *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

When shall Sorbact® Compress not be used?

There are no known contraindications to the use of Sorbact® Compress. Warnings and precautions should be noted.

Warnings and Precautions

Do not use if you have previously experienced hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate).

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of bacteria.

Do not re-use. Sorbact® Compress is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection.

Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed.

Do not wrap around fingers and toes, as it may impair circulation.

Do not re-sterilize.

Consult a healthcare professional if:

- Wound fails to improve or deteriorates
- Side effect is observed
- Signs of infection such as red, warm, swollen or painful skin appear

How to use Sorbact® Compress

1. Clean the wound and dry the surrounding skin.

2. Open the packaging and apply Sorbact® Compress without touching the dressing surface that will be in contact with the wound. The dressing can be unfolded and may overlap the wound margins if needed. Make sure that the product comes in direct contact with the complete wound surface to allow bacteria to bind to the dressing.

3. If the dressing is cut, use clean scissors. Do not tear the dressing. Discard any unused dressing.

4. Apply a secondary dressing appropriate for the exudate level such as an absorbent pad or gauze, and fixate.

5. The dressing change frequency depends on the condition of the wound and the exudate level. The dressing should not be left in place for more than 7 days.

6.Should the dressing adhere to the wound, moisten the dressing with water or saline solution to assist removal and to avoid disruption of the healing wound.

Storage and Disposal

Sorbact® Compress shall be stored dry and kept away from sunlight.

Disposal should be made according to local environmental procedures.

Notice to User

Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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ABIGO Medical AB, Vapenvägen 1, SE-696 33 Askersund, Sweden