
Sorbact® is a registered trademark of ABIGO Medical AB


Date of issue: 2022-05  
Item code: JBR18089.01 / 16571-90109-00

EN Instructions for Use  
SV Bruksanvisning  
FI Käyttöohje  
DA Brugsanvisning  
NO Brukerveiledning  
FR Mode d'emploi

ET Kasutusjuhised  
LV Lietošanas norādījumi  
LT Naudojimo instrukcijos  
PL Instrukcja użyczkowania  
TR Kullanım Talimatları  
FA سوئیچلیل‌های نحوه استفاده

ABIGO Medical AB, Väpenvägen 1, SE-696 33 Askersund, Sweden  
+ 46 31 748 49 50 • sorbact.com

## EN

Read the instructions for use before starting treatment with Sorbact® Surgical Dressing.

### Device Description

Sorbact® Surgical Dressing is a bacteria and fungi binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACC™-coated wound contact layer placed on an absorbent pad and a transparent backing film with an acrylic adhesive. The shower- and bacteria-proof backing film protects against external contamination and allows excess fluid to evaporate. Sorbact® Surgical Dressing prevents and treats wound infections and facilitates the wound healing process. Sorbact® wound contact layer is shown to bind wound related pathogens up to seven days *in vitro*.

### Intended Purpose

Sorbact® Surgical Dressing is intended for use in management of clean, contaminated, colonized or infected dry to low exuding wounds, such as surgical wounds and traumatic wounds.

### Intended User and Use Environment

Sorbact® Surgical Dressing is intended to be used on children and adults. Sorbact® Surgical Dressing is intended to be used by healthcare professionals, or by lay persons under supervision of a healthcare professional. The dressing is intended to be used in healthcare facilities and home environment.

### Mode of Action

Sorbact® Surgical Dressing binds common wound microorganisms, such as *Staphylococcus aureus* (including MRSA), *Streptococcus spp.*, *Escherichia coli*, *Pseudomonas aeruginosa* and *Candida albicans*, as shown *in vitro*. There is no release of antimicrobial agents to the wound.

### Clinical Benefits

Use of Sorbact® Surgical Dressing can reduce bioburden and has been proven to prevent infection. The dressing helps reducing odor and pain, and can improve wound healing and decrease wound size.

### Contraindications

There are no known contraindications to the use of Sorbact® Surgical Dressing. Warnings and precautions should be noted.

### Warnings and Precautions

Do not re-use. Sorbact® Surgical Dressing is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection.

Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed.

Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACC™], cellulose acetate, polyurethane, polyacrylate, polypropylene and viscose).

In common with adhesive products, cases of irritation to the surrounding skin may occur. It should also be noted that inappropriate use or too frequent dressing changes, particularly in patients with fragile skin, may result in superficial skin trauma.

Do not use in combination with fatty products, such as ointments, creams and solutions as they may decrease the binding of microorganisms.

Do not re-sterilize.

If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

### Instructions for Use

1. Prepare the wound and surrounding skin according to local clinical practice. Ensure that the surrounding skin is clean and dry.

2. Select an appropriate dressing size for the wound. Do not cut the dressing.

3. Remove the dressing from the peel pouch using an aseptic technique.

4. Remove the adhesive protective liners marked (1) from the wound contact side of the dressing and apply the dressing over the wound. Ensure that the Sorbact® wound contact layer comes into direct contact with the complete wound surface.

5. Press the borders of the dressing to the surrounding skin and remove the stabilizing film marked (2) on top. Start at the green tab.

6. The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.

7. For easy removal, it is recommended to remove the dressing by pulling and stretching the film of the dressing parallel to the skin.

Should the dressing adhere to the wound, moisten the dressing to assist removal and to avoid disruption of the healing wound. Sorbact® Surgical Dressing is MRI safe.

Prior to radiation therapy, remove Sorbact® Surgical Dressing if placed in the radiation field. A new dressing can be applied following treatment.

**Storage and Disposal**  
Sorbact® Surgical Dressing shall be stored dry and kept away from sunlight.

Disposal should be made according to local environmental procedures.

**Notice to User**  
Any serious incident that has occurred in relation to the device should be reported to ABIGO Medical AB (complaints@abigo.com) and the competent authority of your state.

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