3M™ Littmann® Snap Tight Soft-Sealing Eartips Commonly Asked Questions

Question:

Which type of eartips do I need?

Answer:

Ribbed eartubes on all Littmann stethoscopes manufactured after May 1994 offer a "Snap-tight" fit that will not loosen during routine use. Littmann stethoscopes manufactured before May 1994 feature a threaded eartip design.

Please note: "Threaded" eartips have been discontinued and are no longer available.

Question:

How do I determine the color of eartips to order:

Answer:

Black is used for black tubed stethoscopes. Gray is used for gray or colored tubed stethoscopes. Firm tips are available in gray only.

Question:

What is the charge for parts for my 3M[™] Littmann® stethoscope?

Answer:

Please call the 3M Health Care Service Center at 1-800-292-6298 for item pricing. The Service Center accepts Visa, Mastercard and American Express. All orders are shipped by First Class Mail.

For More Information

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



Health Care