

MULTIGATE

Interpose®



Non-Adherent Absorbent Dressings

Interpose® & Interpose® Lite Range

Our Interpose® dressings are ideal for open wounds that have light to moderate exudate. The dressing consists of a central absorbent layer made of natural cotton which is heat laminated on both sides by a net-like film. The high density polyethylene film provides controlled and consistent porosity which helps to draw exudate from the wound. The dressings' non-adherent properties help to prevent it from sticking to the wound as it heals, resulting in little to no pain being experienced by the patient during removal.



- ✓ Protects and promotes healing
- ✓ Non-adherent to preserve integrity of delicate skin on removal
- ✓ Convenient peel pack (individually wrapped)
- ✓ Australian TGA listing AUST L 28728

Non-Adherent Absorbent Dressing (Peel Pack of 1)

RE-ORDER	DESCRIPTION	SIZE	STERILE / NON-STERILE	PACK PER BOX	BOX PER CARTON	PACK PER CARTON
INTERPOSE®	(Peel Pack of 1)					
30-381	Interpose® Non-Adherent Absorbent Dressing	5cm x 5cm	Sterile	100	6	600pks
30-382	Interpose® Non-Adherent Absorbent Dressing	10cm x 10cm	Sterile	100	4	400pks
30-383	Interpose® Non-Adherent Absorbent Dressing	10cm x 20cm	Sterile	100	2	200pks
30-385	Interpose® Non-Adherent Absorbent Dressing	25cm x 25cm	Sterile	-	-	120pks
30-389	Interpose® Non-Adherent Absorbent Dressing	45cm x 90cm	Sterile	-	-	25pks
INTERPOSE® Lite						
30-501	Interpose® Lite Non-Adherent Absorbent Dressing	7.5cm x 5cm	Sterile	100	8	800pks
30-502	Interpose® Lite Non-Adherent Absorbent Dressing	7.5cm x 7.5cm	Sterile	100	6	600pks
30-503	Interpose® Lite Non-Adherent Absorbent Dressing	7.5cm x 10cm	Sterile	100	6	600pks
30-504	Interpose® Lite Non-Adherent Absorbent Dressing	7.5cm x 20cm	Sterile	100	4	400pks

For more information please contact your local Business Development Manager:

(02) 8718 2888 or visit multigate.com.au