

## Universal Securement Device

Our Universal Securement Device is designed to optimise the securement of indwelling catheters while ensuring patient comfort. An adhesive strip fastener and soft material tabs provide two levels of securement to ensure the catheter stays in place.

- ⊗ Maximise patient comfort
- ⊗ Two levels of securement to avoid sliding and dislodgement
- ⊗ Low irritant adhesive dressing is soft and comfortable on the skin
- ⊗ Re-adjustable soft material tabs
- ⊗ Simple to apply and remove



### Universal Securement Devices

STERILE

RE-ORDER	DESCRIPTION	DIMENSION	SIZE	CARTON QTY
UR060000	Universal Securement Device	6.5cm x 10.5cm	Standard	50
UR060001	Universal Securement Device	8.5cm x 12.5cm	Large	50